

## **Legal Notice**

© copyright 2015 Romaco Innojet GmbH

This complete presentation is property of Romaco Innojet GmbH in D-Steinen.  
Publishing in any media neither of the total document nor of parts of it  
like pictures or supply to third parties  
are not permitted without written authorization.

Diese komplette Präsentation ist Eigentum der Romaco Innojet GmbH in D-Steinen.  
Sowohl die Veröffentlichung in Medien als auch die Weitergabe,  
weder des gesamten Dokuments noch von Auszügen davon,  
wie Bildern, sind nicht gestattet, ohne vorherige, schriftliche Genehmigung.

# PHARMA*process*

Innovation Forum in Pharmaceutical Process

**High Sophisticated Technologies  
for Oral Solid Dosage Forms:**

**Hotmelt Granulation and Coating**

**Part 1 – by Romaco Innojet**

**Latest Fluid Bed Systems for Granulation and  
Coating for Hotmelt Applications**

**Part 2 – by Hermes Pharma**

**R&D of High Sophisticated OSD and  
Hot Melt Coating Processing**

[www.pharmaprocessforum.com](http://www.pharmaprocessforum.com)

## Part 1

### Latest Fluid Bed Systems for Granulation and Coating for Hotmelt Applications



**Quality By Design approach is essential:**

**To chose the right technology and the corresponding equipment  
for R&D from the Laboratory to  
Pilot Plant – Production Scale**

© 2015 by Romaco Innojet GmbH



## Hot-Melt-Coating (HMC)

**which is: Spraying with melted materials (Granulation, Coating)**

- fast up-coming as a high profitable solution for production in Food and Pharma
- looking to avoid the shortfalls of solvent-based granulation and coating
- to better control taste (masking) and release (efficacy) parameters
- requires special know-how and insight expertise to develop the suitable technical solution

## Obstacles


- there are pretty rare so called Ready-To-Use materials
- just few or no materials registered for Pharma and Food applications in comparison to Polymer Coatings (aqueous or for organic solvents)
- conventional Fluid Bed Devices (Top-Spray or Bottom Spray) are very limited or even unable to perform a robust, reproducible and high efficient granulation or coating process due to their physical characteristics

## Hot-Melt-Coating (HMC)

### Functionality Of HMC-System

#### 3 Major Functional Core Components to achieve required physical conditions

1. **Fluid Bed** for high precisely controlled and homogeneous particle movement (int. pat Romaco Innojet ORBITER)
2. **Spray Nozzle** (int. pat. Romaco Innojet ROTOJET) – controlled, high precisely spray supply
3. **In-Process Filters** (int. pat. Romaco Innojet SEPAJET) – permanent product recovery
4. Melting Device with heated pump and dosage system





## Hotmelt Advantages: Process Time in a Romaco Innojet V 400 (400 ltr. process batch)

Process data	Hot Melt process	Standard process
Process time	1h	5 h
Quantity of air	6000 m <sup>3</sup> /h	6000 m <sup>3</sup> /h
Inlet air temp.	Not heated	70 °C
Melt temp./ liquid temp.	90 °C	20 ° C
Temp. of spray air	90 °C	60 °C

Calculation basis:

$c_p$  wax/fat: 2,7 J/g/K

melt enthalpy H wax/fat: 160 J/g



## Hotmelt Advantages: Energy Consumption in a Romaco Innojet V 400 (400 ltr. process batch)

	Hot Melt process	Standard process
Heating up the inlet air	- - -	5 h x 125 kW = 500 kWh
Heating up the melt/wax	10,5 kWh	- - -
Melting energy wax	200 kg wax = 8,8 kWh	- - -
Spray air heating	1 h x 15 kW = 15 kW	5 h x 10 kW = 50 kWh
Tube/ house/ nozzle heating	1 h x 2 kW = 2 kW	- - -
<b>Total</b>	<b>36, 3 kWh</b>	<b>550 kWh</b>



## Hotmelt: Advantages

- No liquid / moisture in the system
- No evaporation of aqueous solution
- Very short process times
- Green Pharma – no organic solvents – natural materials

### Application example:

HMC – Taste masking of granules

Starter: 250 g granules

50 % coating -> 125 g fat/wax

Spray rate: 7 g/min

**Process time: 18 min**

Taste masking with polymer dispersion:

15% polymer in dispersion

125 g polymer – 833 g dispersion

Spray rate: 7 g/min


**Process time: 120 min**

**(5 g/min - 170 min)**






## Product Sample: Taste masking of bitter tasting API



Pin Head



Process parameters

Starter material:  
500 g crystals



Final product :  
1000 g coated crystals

Process time: 115 min  
Spray rate: 3,6 – 8 g/min



## Product Sample: Malic Acid Crystals protected against humidity

- +30% fat, size 200 $\mu$ m – 500 $\mu$ m



## Product Sample: Encapsulation of Alginate Capsules made by Coacervate Method

### Encapsulation: 1<sup>st</sup> pre-drying in Ventilur<sup>®</sup> System – 2<sup>nd</sup> Hotmelt Coating

- liquid API-s
- flavors
- to provide handling of liquids as a dry product





## Application Example: Encapsulation of Probiotics

### Unique Multiple-Layered Micro-Encapsulation – a Single-Pot-Process

- started with granulation and completed with multi-layer coating – non-stop just the spray nozzle got changed while product was under control
- to provide Pro-Biotic Bacteria with maximum heat resistance
- to ensure highest biological efficacy in the lower GIT
- to achieve Healthy Food Comprising THEM 25847-WO-09

**Yellow** – probiotics granules (the core)

**Blue** – moisture barrier coating layer

**Pale blue** – subcoating-intermediate layer

**Brown** – acid-resistant layer (enteric coating)

**Pale brown** – heat-protection layer





## Romaco Innojet VENTILUS® 3-in-1 Drying, Granulation & Coating System

### Laboratory scale

- V 2.5 - 1
- V 5





### Pilot scale

- V 10
- V 25
- V 50




### Production scale

- V 600
- V 800
- V 1000



## Romaco Innojet VENTILUS® Granulation & Coating System







## Romaco Innojet VENTILUS® 3-in-1 System for Drying, Granulation & Coating



int. pat., developed by Dr. h.c. Herbert Hüttlin



## Romaco Innojet VENTILUS® Hotmelt device IHD



Laboratory scale

VENTILUS® equipment is working with only one centrally placed spray nozzle ROTOJET Hotmelt which forms the basis for successful Hotmelt Applications




Production scale



## Romaco Innojet VENTILUS®

1<sup>st</sup> Controlled Particle Movement – Homogeneous Temperature Supply  
Booster ORBITER type IBO




Example: Production scale unit 800 liters batch size

int. pat., developed by Dr. h.c. Herbert Hüttlin



## Romaco Innojet VENTILUS®

1<sup>st</sup> Controlled Particle Movement – Homogeneous Temperature Supply  
Booster ORBITER type IBO




Romaco Innojet booster ORBITER (side view, open air gaps highlighted).



air guiding blades – rear side

## Romaco Innojet VENTILUS®

1<sup>st</sup> Controlled Particle Movement – Homogeneous Temperature Supply  
Booster ORBITER type IBO



The relative number of air-gaps per Booster surface is constant


- ▶ The relative air-quantity per product-volume is constant
- ▶ Nearly linear up-scaling characteristic

int. pat., developed by Dr. h.c. Herbert Hüttlin




## Romaco Innojet VENTILUS®

1<sup>st</sup> Controlled Particle Movement – Homogeneous Temperature Supply  
Booster ORBITER type IBO



Romaco Innojet booster ORBITER:  
Horizontal – vertical air stream





Sieve, Conidur and slot bottoms:  
Only vertical or tangential air stream

The prolonged path of the process air on the booster  
ORBITER enhances **the drying efficiency by 25 %**

## Romaco Innojet VENTILUS® Functional principle of booster ORBITER

The product movement follows a **curve** without any impact on the container wall

The curve is created by “**air guiding pins**” below the booster plates






## Romaco Innojet VENTILUS® Booster ORBITER



## Romaco Innojet VENTILUS® 2<sup>nd</sup> Controlled and Precise Spray Supply Spray nozzle ROTOJET



Detailed view on spray liquid gaps



Individual ROTOJET parts



Laboratory- and  
pilot scale  
(e.g. type INR  
10/25)



Production scale  
(type INR  
50/75/100/125)

int. pat., developed by Dr. h.c. Herbert Hüttlin




## Romaco Innojet VENTILUS®

### 2<sup>nd</sup> Controlled and Precise Spray Supply

### Comparison spray nozzle ROTOJET vs. conventional nozzles


ROTOJET (all sizes)

Two-substance spray nozzles





## Romaco Innojet VENTILUS® 2<sup>nd</sup> Controlled and Precise Spray Supply Scale-Up of spray nozzle ROTOJET




Spray gap geometry remains the same, independent of the nozzle size

► Scale-up is the basic philosophy

int. pat., developed by Dr. h.c. Herbert Hüttlin


## Romaco Innojet VENTILUS® Spray nozzle ROTOJET





## Romaco Innojet VENTILUS®

Functional principle of booster ORBITER and spray nozzle ROTOJET




**Product movement** and **nozzle spray** work into the same direction

int. pat., developed by Dr. h.c. Herbert Hüttlin



## Romaco Innojet VENTILUS®

Functional principle of booster ORBITER and spray nozzle ROTOJET



- Spraying and moving **from the center** make the particles disperse
- Very high spray rates are possible

## Romaco Innojet VENTILUS®

### 3<sup>rd</sup> Controlled and Permanent Powder Recovery System SEPAJET



view into the filter dome from working position with product container opened



simple and fast change of filter set

int. pat., developed by Dr. h.c. Herbert Hüttlin






## Romaco Innojet VENTILUS®

### 3<sup>rd</sup> Controlled and Permanent Powder Recovery System SEPAJET




Stainless steel carrier basket  
with textile filters





int. pat., developed by Dr. h.c. Herbert Hüttlin

## Romaco Innojet VENTILUS®

### 3<sup>rd</sup> Controlled and Permanent Powder Recovery System SEPAJET Comparison conventional filters vs. SEPAJET



SEPAJET filter with  
cylindrical product container



Conventional filter with conical  
product container



- Conventional filters accelerate the process air at the filter tubes' lower edge or in single chamber systems the process is stopped for filter cleaning/shaking
  - Fine material stays in the filter
  - Non-stop and permanent powder recovery and filter cleaning



## Romaco Innojet VENTILUS® Powder recovery system SEPAJET






## Romaco Innojet VENTILUS®

### 3<sup>rd</sup> Controlled and Permanent Powder Recovery System SEPAJET

#### Comparison conventional filters vs. SEPAJET


Inner surface area in contact with the product (including filters)

Conventional Double Chamber Tube-Filter system



Inner surface area: FBD 60: 50 m<sup>2</sup>

SEPAJET Systems



V 50: 8 m<sup>2</sup>

V 100: 12 m<sup>2</sup>

## Romaco Innojet VENTILUS®

### Laboratory System – Pilot Plant – Production Scale

Nearly linear up-scaling from laboratory up to production scale

#### **Excellence due to Physical Features of Patented Design**

- + 99% dust free granules
- + homogeneous film coating
- + exactly controlled process
- + hot melt coating
- + reduced process time
- + increased yield





## Romaco Innojet VENTILUS® Laboratory System – Pilot Plant – Production Scale





## Romaco Innojet VENTILUS®

Laboratory System – Pilot Plant – Production Scale



Machine opened

GMP compliant access to inner parts of the machine

Contamination free loading – unloading

Discharge container (closed discharge)



## Romaco Innojet VENTILUS® Laboratory System – Pilot Plant – Production Scale



GMP compliant access to inner parts of the machine



GMP compliant change of spray nozzle – available in-process



## Romaco Innojet VENTILUS® Laboratory System – Pilot Plant – Production Scale




## Romaco Innojet VENTILUS®

Laboratory System – Pilot Plant – Production Scale

### 3-in-1 Multi Purpose System


1. Granulation – 2. Film Coating Micro Tablets – 3. Pellets Layering and Coating



Endproduct



Starterpellets



Micro Tablets







## Romaco Innojet VENTILUS®

Laboratory System – Pilot Plant – Production Scale

### EXAMPLE:


**MUPS – Multi Unit Particles Systems implemented in RDT or in Capsules, Stick Packs**



Capsules



Stickpacks/  
Sachets



Micro CT of MUPS tablets

© copyright & all rights reserved  
by ACINO PHARMA AG 2010 in CH-Basel



## Romaco Innojet VENTILUS®

Laboratory System – Pilot Plant – Production Scale

### EXAMPLE: Film Coating, Crop Protection

Flower Seeds Coating

- homogeneous & fast
- gentle & reproducible
- no agglomerations




## Romaco Innojet VENTILUS®

Laboratory System – Pilot Plant – Production Scale

### EXAMPLE: Solidification and Unique Encapsulation and Preservation Multiple-Layered Micro-Encapsulation

- to provide sensitive natural oils with a superior stability
- preservation against oxygen and humidity
- improved stability for extended shelf life




Natural Liquid Oil

Oxygen and Humidity Barrier

Seal Coating Layer

Intermediate Coating Layer



- Flavoring Agents (e.g. Lemon Oil, Peach Oil, etc.)
- Poly-Unsaturated Fatty Acids (PUFA , e.g. Omega 3, etc.)

© 2012 by Dr. Adel Penhasi SPAI/DeGama Group

## **Romaco Innojet VENTILUS®**

Granulation & Coating – Aqueous – Organic Solvents – Hotmelt Applications

**Fast and Reproducible – Direct Granulate  
Single-Pot Processing, yield >98% of batch**

1<sup>st</sup> Granulation

2<sup>nd</sup> in-process change of the single Spray Nozzle

3<sup>rd</sup> application of a sub-coating

4<sup>th</sup> application of Enteric Coating

5<sup>th</sup> after total process time: Final sieving of  
ready-made product for packing into Stick Packs





## Contact:

**ROMACO INNOJET**  
*beyond technology* *Herbert Hüttlin*

e-mail: [innojet@romaco.com](mailto:innojet@romaco.com)

Romaco Innojet GmbH  
Daimlerstr. 7  
79585 Steinen | Germany

Int. Pat. Granulation & Coating Technology  
[www.romaco.com](http://www.romaco.com)





Get the dose right™

## Legal Notice

© copyright 2015 Hermes Arzneimittel GmbH – Division HERMES PHARMA

© 2015 by HERMES PHARMA Dr. Detlev Haack

This complete presentation is property of Hermes Arzneimittel GmbH,  
Division HERMES PHARMA in D-Pullach (München).

Publishing in any media neither of the total document nor of parts of it  
like pictures or supply to third parties  
are not permitted without written authorization.

Diese komplette Präsentation ist Eigentum der Hermes Arzneimittel GmbH,  
Division HERMES PHARMA in D-Pullach (München).

Sowohl die Veröffentlichung in Medien als auch die Weitergabe,  
weder des gesamten Dokuments noch von Auszügen davon,  
wie Bildern, sind nicht gestattet, ohne vorherige, schriftliche Genehmigung.

## Part 2

### R&D of High Sophisticated OSD and Hot Melt Coating Processing

**HERMES**  
PHARMA

*Get the dose right™*

### New Taste Masking Technology to Facilitate User-Friendly Pharmaceuticals

Hermes Pharma, INNOJET Herbert Hüttlin,  
RCPE and Karl-Franzens University

© 2015 by HERMES PHARMA Dr. Detlev Haack



## Aim

## Aim

- Development of taste masked powder or granulate with an immediate release profile
- Packaging and sales in a Stick Pack

**HERMES  
PHARMA**

Get the dose right™

## Why Oral Solid Dosage in a Stick Pack?

- Medicine-To-Go: Convenient Stick Pack
- Ready and easy to use – taking without water
- Faster liberation than tablets:  
Improved bio-availability
- Easy to swallow – pediatrics, geriatrics






## Hurdles and Approach

**HERMES**  
PHARMA

Get the dose right™

### Hurdles

- Amount of single dose approx. 2 g (salivation sufficient)
- Taste masking with matrix not enough
- Amount of matrix limited



### Approach

- Hot Melt Coating of API-crystals with lipid excipients  
→ taste masking with immediate release profile





## Hot Melt Coating with Lipid Excipients

**HERMES  
PHARMA**

Get the dose right™

### Advantages

- Short process times (saving of solvent evaporation step)
- Coating materials toxicologically uncritical and from natural sources (Green Pharmacy)
- Dissolution pH-independent (compatible with acidic matrix components)
- Lipids: Pleasant taste, odor and feeling in the mouth
- Applicable for processing with implemented PAT

### Disadvantages

- Polymorphic changes during storage
- Retardation of dissolution profile (sustained release)



Romaco Innojet VENTILUS® Laboratory System V 2.5/1

© int. pat. by Dr. h.c. Herbert Hüttlin



# QTPP Formulation Development Strategy

## Quality Target Product Profile

**Storage Stability:** Climate zone IV b  
(30° C, 75 % r.h.) 12 months  
**Fast Dissolution:** 85 % in 30 min  
**Taste Masking:** 1 min

## Storage Stability:

Avoiding of polymorphic changes

→ Adding of modifiers (emulsifiers):

- Induction of the stable  $\beta$ -form:

Increasing of emulsifier content 10 % to 30 %

(21 h vs. 2 h to achieve stable  $\beta$ -form)

→ Additional: lower process temperatures applicable



## DoE

**HERMES**  
PHARMA

Get the dose right.™

## Screening DoE (Frac Fac Res IV) for linear interactions

## Definition of Critical Process Parameters as Input Parameters:

### Input parameter:

- Spray Rate (SR)
- Spray Pressure (SP)
- Coating Amount (CA)
- Emulsifier Content (EMU)
- Air Flow Rate (Air)



## DoE


### Interaction with Critical Quality Attributes

**HERMES**  
PHARMA

Get the dose right™

#### Output parameter: Dissolution rate

- **Emulsifier Content (EMU):**  
Increasing EMU:  
faster dissolution  
and shorter transformation  
time to the stable  $\beta$ -form
- **Coating Amount (CA):**  
Increasing CA:  
Slower dissolution



R<sup>2</sup>: 0.87  
Q<sup>2</sup>: 0.82  
Model validity: 0.57  
Reproducibility: 0.99





## DoE


### Interaction with Critical Quality Attributes

**HERMES**  
PHARMA

Get the dose right™

**Output parameter:** Taste masking efficiency (in vitro pH measurement)

- **Coating Amount (CA)**
- **Air flow rate – quadratic term ( $\text{Air}^2$ )**
- **Spray Pressure (SP)**  
Increasing CA/Air/SP:  
improved taste masking
- **Emulsifier Content (EMU)**  
Increasing EMU:  
degraded taste masking



R<sup>2</sup>: 0.92


Q<sup>2</sup>: 0.86

Model validity: 0.95

Reproducibility: 0.81

HERMES  
PHARMA

Get the dose right™

[illegible]

SEM picture of API-Crystal coated with a lipid/emulsifier mixture

54 © 2015 by HERMES PHARMA Dr. Detlev Haack




## Stability Study

### Dissolution immediate release

Specification Limits:  
85% in 30 min

**HERMES**  
PHARMA

Get the dose right™





## HMC – Summary – Conclusion

**Hot-Melt Coating: An available and proved alternative technology to replace conventional Polymer Coatings**

**HERMES  
PHARMA**

Get the dose right™

### Process Technology

- Very short processing times in comparison to conventional coatings
- Robust and reproducible process
- Easy/simple scale-up
- PAT solutions to control and determine the product quality
- No window fouling on probe in pilot and production plant
- PAT controls the endpoint of the coating process in real time





## HMC – Summary – Conclusion

**Hot-Melt Coating: An available and proved alternative technology to replace conventional Polymer Coatings**

**HERMES  
PHARMA**

Get the dose right™

### Pharmaceutical Aspects

- Taste masking for user friendly **Oral Solid Dosage** forms achievable
- Excellent bio-availability
- Increased customers convenience
- Stable product over shelf life
- Immediate release profile and reliable taste masking
- Significant cut of business costs and increase of profitability
- Solvent free – Green Pharma
- Sustained and controlled release (Enteric Coating) achievable

## Contact:

Dr. Detlev Haack – Director R&D

e-mail: [haack@hermes-pharma.com](mailto:haack@hermes-pharma.com)

Hermes Arzneimittel GmbH  
Division HERMES PHARMA  
Georg-Kalb-Strasse 5-8  
82049 Pullach | Germany

**HERMES**  
PHARMA

*Get the dose right™*

Full Service R&D – Contract Manufacturing  
[www.hermes-pharma.com](http://www.hermes-pharma.com)